

# *Applications of the Crystallographic Search and Analysis System CRYSTDAT in Materials Science*

---

Volume 94

Number 1

January-February 1989

---

## **T. Siegrist**

AT&T Bell Laboratories  
600 Mountain Avenue  
Murray Hill, NJ 07974

Numerical database systems have recently become available online. Their enhanced search capabilities and fast retrieval of data make them a valuable tool in research. In particular, CRYSTDAT which is a search and analysis system for NBS CRYSTAL DATA has proven to be powerful in the identification of crystalline materials. In con-

junction with a single-crystal x-ray diffractometer, a qualitative as well as quantitative phase determination is easily performed. The use of CRYSTDAT will be illustrated in several examples.

**Key words:** crystallographic database; CRYSTDAT; high  $T_c$  superconductors.

---

## **Introduction**

Numerical database systems compiling inorganic crystallographic data are well established. Among the best known are the JCPDS Powder Diffraction File [1], Structure Reports (Strukturbericht) [2] and Crystal Data Determinative Tables [3]. To retrieve data from the databases in printed form, elaborate and sometimes limited and inconvenient search methods are required. In contrast, online computer systems are highly efficient. CRYSTDAT [4], which is a search and retrieval system for NBS CRYSTAL DATA [5] is accessible online and makes use of efficient search algorithms. It is therefore possible to search many parameter fields efficiently and the retrieved data can be quickly evaluated. In the following, several examples on the possible use of CRYSTDAT will be given. These examples only cover a small part of the different search routines available and use a subset of the options implemented.

## **Applications**

Currently, NBS CRYSTAL DATA contains unit cell, composition, formula, symmetry, and reduced cell data of approximately 115,000 inorganic and organic compounds up to 1985. Using CRYSTDAT, each of these data fields can be searched individually or they can be combined through logical operators. The examples will be restricted to inorganic oxides, illustrating possible uses of CRYSTDAT in the research on high  $T_c$  ceramic superconductors. In particular, the first example will show how the database was used in the study of the Ba-Y-Cu-O system. After logging into the account and loading the database, the system replies:

\_\_\_\_\_

Copyright U.S. Department of Commerce,

National Bureau of Standards, 1985

---

The coverage is for the years 1900–1985

```
>find class i (find inorganic compounds)
--Set 1 created with 59951 hits
>set limits 1 (limits following searches to set 1)
--Limit Set=1
>find ele O (find oxides)
--Set 2 created with 31062 hits
>set limits 2
--Limit Set=2 (limits following searches to set 2)
```

```
>find ele Ba.and.Cu.and.3 (search for ternary barium cuprates)
--Set 3 created with 3 hits
>show 3 (present the results of set 3)
```

- 1

-- 2

<sup>1</sup> Output abbreviations are as follows: ID=compound identification no., RC=reduced cell, CD=Crystal Data, EM=empirical formula, FO=chemical formula, NM=name, AC=author's cell, RF=journal reference.

ID : 800744  
 RC : a=15.82 b=15.82 c=15.82 al=109.5 be=109.5 ga=109.5  
 CD : sys=cubic spgr(CD)=Im3m spno=229 den=5.6(g/cc) z=90  
 EM: Ba Cu O2  
 FO : Ba Cu O2  
 NM: Barium dioxocuprate  
 AC : a=18.27 spgr(A)=Im3m  
 RF : Z. Naturforsch. B,32,121,1977

-- 3

#### b. Search for Ternary Cu-Y-Oxides

>find Y.and.Cu.and.3  
 --set 4 created with 2 hits  
 >show 4

ID : 711441  
 RC : a=3.50 b=10.80 c=12.46 al=90.0 be=90.0 ga=90.0  
 CD : sys=orthorhombic spgr(CD)=Pn21a spno=33 den=5.4(g/cc) z=4  
 EM: Cu2 O5 Y2  
 FO : Cu2 Y2 O5  
 NM: Copper Yttrium Oxide  
 AC : a=10.799 b=3.4960 c=12.456 spgr(A)=Pna21  
 RF : 00GRNT,,,1981

-- 1

ID : 811630  
 RC : a=3.52 b=3.52 c=11.42 al=90.0 be=90.0 ga=120.0  
 CD : sys=hexagonal spgr(CD)=P63/mmc spno=194 den=5.0(g/cc) z=2  
 EM: Cu O2 Y  
 FO : Cu Y O2  
 NM: Copper(i) yttrium oxide  
 AC : a=3.5206 c=11.418 spgr(A)=P63/mmc  
 RF : J. Solid State Chem.,49,232,1983

-- 2

#### c. Search for Ternary Ba-Y-Oxides

>find ele Y.and.Ba.and.3  
 --Set 5 created with 4 hits  
 >show 5

ID : 705238  
 RC : a=3.45 b=10.39 c=12.11 al=90.0 be=90.0 ga=90.0  
 CD : sys=orthorhombic spgr(CD)=Pnab spno=60 den=7.6(g/cc) z=0  
 EM: Ba O4 Y2  
 FO : Ba Y2 O4  
 NM: Barium Yttrium Oxide  
 AC : a=10.388 b=12.110 c=3.448 spgr(A)=Pnab  
 RF : Mater. Res. Bull.,9,1631,1974

-- 1

ID : N108141 -- 2

RC : a=3.45 b=10.41 c=12.12 al=90.0 be=90.0 ga=90.0

CD : sys=orthorhombic spgr(CD)=Pnam spno=62 den=5.8(g/cc) z=4

EM: Ba O4 Y2

FO : Ba Y2 O4

NM: Barium yttrium oxide (1<sup>-2</sup> 2<sup>-4</sup>)

AC : a=10.415 b=12.120 c=3.455 spgr(A)=Pnam

RF : Z. Naturforsch. B,19,955,1964

ID : 705239 -- 3

RC : a=4.38 b=4.38 c=11.85 al=90.0 be=90.0 ga=90.0

CD : sys=tetragonal den=8.2(g/cc) z=0

EM: Ba2 O5 Y2

FO : Ba2 Y2 O5

NM: Barium Yttrium Oxide

AC : a=4.3771 c=11.852

RF : Mater. Res. Bull.,9,1631,1974

ID : 706037 -- 4

RC : a=6.11 b=6.11 c=25.17 al=90.0 be=90.0 ga=120.0

CD : sys=hexagonal den=2.0(g/cc) z=0

EM: Ba3 O9 Y4

FO : Ba3 Y4 O9

NM: Barium Yttrium Oxide

AC : a=6.1102 c=25.172

RF : Mater. Res. Bull.,9,1631,1974

#### d. And Finally Search for the Quaternary Ba-Y-Cu-Oxides

>find ele Ba.and.Y.and Cu.and.4

--Set 6 created with 1 hits

>show 6

ID : 809291 -- 1

RC : a=5.66 b=7.13 c=12.18 al=90.0 be=90.0 ga=90.0

CD : sys=orthorhombic spgr(CD)=Pbnm spno=62 den=6.2(g/cc) z=4

EM: Ba Cu O5 Y2

FO : Y2 Ba Cu O5

NM: Dyttrium barium copper oxide

AC : a=7.132 b=12.181 c=5.658 spgr(A)=Pbnm

RF : J. Solid State Chem.,43,73,1982

Now we can (attempt to) construct a preliminary phase diagram of the Ba-Y-Cu-O system that serves as the starting point of the analysis of this phase space.

#### e. Identification of a Crystalline Material

In the course of studying the phase system, small black crystals were obtained. From 6 reflections measured on the CAD4 diffractometer, a unit cell of  $3.495 \times 6.23 \times 10.795 \text{ \AA}$  was inferred. This unit cell serves as the starting point for a database search to identify the phase. In this case, the reduced cell is searched with a preset tolerance to retrieve all phases with similar lattice parameters.

>cells (calculate the reduced cell and set up the search parameters)

---Input Cell---

a ? > 3.495

b ? > 6.23 (if blank, default is the previous value)

c ? > 10.795

alpha ? > (the default is 90 °)

beta ? > (if blank, default is the previous value)

gamma ? >

Lattice Type ? (default P) >

Tolerance ? (default 0.10) >

```
-----
>Input cell+vol.      3.50   6.23   10.80   90.00   90.00   90.00   235.05
>Lattice type         P
>Red. Cell+Vol        3.49   6.23   10.79   90.00   90.00   90.00   235.05
>Niggli Matrix        12.215  38.813  116.532
                      0.000   0.000   0.000
>Metric Symmetry      Orthorhombic
>Tolerance for cell match 0.100
-----
```

Search the data base? (yes/no)>y

--Set 7 created with 0 hits

*The database was searched with the reduced cell of the lattice with edges of 3.49, 6.23, and 10.795 Å. No matches were found. However, the unit cell determined may be too small (a subcell) because some reflections could have been missed on the diffractometer. The database is searched again, but this time, possible supercells are included in the search.*

>cells

---Input Cell---

a ? > 3.495

b ? > 6.23

c ? > 10.795

alpha ? >

beta ? >

gamma ? >

Lattice Type ? (default P) >

Tolerance ? (default 0.10) >

```
-----
>Input cell+vol.      3.50   6.23   10.80   90.00   90.00   90.00   235.05
>Lattice type         P
>Red. Cell+Vol        3.49   6.23   10.79   90.00   90.00   90.00   235.05
>Niggli Matrix        12.215  38.813  116.532
                      0.000   0.000   0.000
>Metric Symmetry      Orthorhombic
>Tolerance for cell match 0.100
-----
```

Search the data base? (yes/no)>n

*The database is not searched with the original cell. Instead, a supercell search [6] will be carried out. CRYSTDAT is instructed to generate possible supercells of 2 × the volume of the input cell. After the supercells have been generated, the database is searched.*

Sub-/Super- cells calculated? (YES/NO)  
 where sub-cell (input cell vol.) / multiplicity  
 and super-cell (input cell vol.) \* multiplicity > y  
 Sub-cell ? (yes/no) > n  
 Super-cell ? (yes/no) > y  
 Multiplicity ? (2/3/4/5/6/7 ... etc) > 2

*Super-cells of 2 times the volume of the input cell*

a	b	c	al	be	ga	vol
3.495	6.230	21.590	90.0	90.0	90.0	470
3.495	12.460	12.464	120.0	90.0	90.0	470
6.230	6.990	11.347	107.9	90.0	90.0	470
6.990	7.143	11.347	81.3	72.1	60.7	470
3.495	10.795	12.460	90.0	90.0	90.0	470
6.990	7.143	10.795	90.0	90.0	119.3	470
6.230	6.990	10.795	90.0	90.0	90.0	470

Search the data base? (yes/no) > y

**\*\*NOW SEARCHING\*\***

*The database is now searched including all the above given supercells of multiplicity 2.*

--Set 8 created with 31 hits

*In this set of retrieved entries, we find the following one:*

ID : 711441 -- 1  
 RC : a=3.50 b=10.80 c=12.46 al=90.0 be=90.0 ga=90.0  
 CD : sys=orthorhombic spgr(CD)=Pn21a spno=33 den=5.4(g/cc) z=4  
 EM: Cu2 O5 Y2  
 FO : Cu2 Y2 O5  
 NM: Copper Yttrium Oxide  
 AC : a=10.799 b=3.4960 c=12.456 spgr(A)=Pna21  
 RF : 00GRNT,,1981

*Supercell 5 from the above list of supercells matched a cell in the database of a compound with the expected elemental composition. In this way, both the correct lattice and an accurate composition of the small black crystals were determined in roughly 30 minutes.*

## 2. Study of the Bi-Sr-Ca-Cu-O System

*In the case of the Bi-Sr-Ca-Cu-O superconductors, the analogous searches are easily done.*

### a. Search for Ternary Bi-Cu Oxides

>find ele Bi.and.Cu.and.3  
 --Set 9 created with 4 hits  
 >show 9

ID : B034138 -- 1  
 RC : a=5.80 b=8.48 c=8.48 al=90.0 be=90.0 ga=90.0  
 CD : sys=tetragonal spgr(CD)=P4/ncc spno=130 den=8.0(g/cc) z=2  
 EM: Bi4 Cu O7  
 FO : Bi4 Cu O7  
 NM: Bismuth copper oxide (4~1~7)  
 AC : a=8.481 c=5.803 spgr(A)=P4/ncc  
 RF : An. Acad. Brasil. Cienc.,38,35,1966

ID : 804574 -- 2  
 RC : a=5.81 b=6.67 c=6.67 al=79.0 be=64.1 ga=64.1  
 CD : sys=tetragonal spgr(CD)=I4 spno=79 den=8.6(g/cc) z=4  
 EM: Bi2 Cu O4  
 FO : Cu Bi2 O4  
 NM: Copper dibismuthate(iii)  
 AC : a=8.484 c=5.813 spgr(A)=I4  
 RF : Z. Anorg. Allg. Chem.,426,1,1976

ID : 704311 -- 3  
 RC : a=5.81 b=8.51 c=8.51 al=90.0 be=90.0 ga=90.0  
 CD : sys=tetragonal spgr(CD)=P4/ncc spno=130 den=8.6(g/cc) z=4  
 EM: Bi2 Cu O4  
 FO : Cu Bi2 O4  
 NM: Copper Bismuth Oxide  
 AC : a=8.510 c=5.814 spgr(A)=P4/ncc  
 RF : C. R. Hebd. Seances Acad. Sci. Ser. C,276,1105,1973

ID : 810281 -- 4  
 RC : a=5.81 b=8.51 c=8.51 al=90.0 be=90.0 ga=90.0  
 CD : sys=tetragonal spgr(CD)=P4/ncc spno=130 den=8.6(g/cc) z=4  
 EM: Bi2 Cu O4  
 FO : Cu Bi2 O4  
 NM: Bismuth copper oxide  
 AC : a=8.510 c=5.814 spgr(A)=P4/ncc  
 RF : Bull. Soc. Fr. Mineral. Cristallogr.,99,193,1976

#### b. Search for Ternary Bi-Sr-Oxides

>find ele Bi.and.Sr.and.3  
 --Set 10 created with 4 hits  
 >show 10

ID : B028466 -- 1  
 RC : a=3.95 b=3.95 c=9.64 al=78.2 be=78.2 ga=60.0  
 CD : sys=rhombohedral den=6.0(g/cc) z=0  
 EM: Bi38 O59 Sr2  
 FO : Bi38 Sr2 O59  
 NM: Bismuth strontium oxide (38~2~59)  
 AC : a=3.952 c=28.09 spgr(A)=R  
 RF : J. Res. Nat. Bur. Stand. Sect. A,68,197,1964

ID : 805401 -- 2  
 RC : a=3.97 b=3.97 c=9.75 al=78.2 be=78.2 ga=60.0  
 CD: sys=rhombohedral spgr(CD)=R-3m spno=166 den=7.8(g/cc) z=9  
 EM: Bi0.76 O1.38 Sr0.23  
 FO : Bi0.765 Sr0.235 O1.383  
 NM: Bismuth strontium oxide (.8/.2/1.1)  
 AC: a=9.75 al=23.49 spgr(A)=R-3m  
 RF : J. Solid State Chem.,35,192,1980

ID : 812013 -- 3  
 RC : a=3.97 b=3.97 c=9.74 al=78.2 be=78.2 ga=60.0  
 CD: sys=rhombohedral spgr(CD)=R-3m spno=166 den=7.8(g/cc) z=4  
 EM: Bi1.72 O3 Sr0.53  
 FO : Bi1.72 Sr0.53 O3  
 NM: Bismuth(iii) strontium(13/4) oxide  
 AC: a=3.971 c=28.41 spgr(A)=R-3m  
 RF : Solid State Ionics,3,457,1981

ID : 709884 -- 4  
 RC : a=4.26 b=9.60 c=9.60 al=87.2 be=77.2 ga=77.2  
 CD: sys=tetragonal spgr(CD)=I4/m spno=87 den=7.2(g/cc) z=0  
 EM: Bi1.10 O2.55 Sr0.90  
 FO : Sr0.9 Bi1.1 O2.55  
 NM: Strontium Bismuth Oxide  
 AC: a=13.239 c=4.257 spgr(A)=I4/m  
 RF : Rev. Chim. Miner.,15,153,1978

#### c. Search for Bi-Cu-Sr or Bi-Cu-Ca Compounds with 4 or 5 Elements

>find ele Bi.and.Cu.and.Sr.and.5  
 --Set 11 created with 0 hits  
 >find ele Bi.and.Cu.and.Sr.and.4  
 --Set 12 created with 0 hits  
 >find ele Bi.and.Ca.and.Cu.and.4  
 --Set 13 created with 0 hits  
 >find ele Bi.and.Cu.and.Ca.and.5  
 --Set 14 created with 0 hits

*No entries containing 4 and 5 elements were found.*

#### d. Find materials with lattices related to that of a superconducting crystal

*From the x-ray study of a superconducting crystal, we obtained a cell of approximately  $5.44 \times 5.43 \times 30.8$  Å. This cell data can be used to retrieve compounds that may be structurally related to the superconductors. For this purpose, tetragonal or orthorhombic cells are retrieved and searched for their c-axis values given by the authors.*

*Carry out several searches to generate SET 17 containing inorganic materials with Bi and O; restrict next search to this set.*

>Find sys O.or.T (orthorhombic or tetragonal systems only)  
 --Set 18 created with 244 hits  
 >Set limits 18 (limits now set to inorganic bismuth oxides belonging to the orthorhombic or tetragonal systems.)



--Limit Set=18

>find cc 28.to.34 (*find compounds with the author's c-axis between 28 and 34 Å*)

--Set 19 created with 8 hits

*In this way, compounds with unit cell c-axes between 28 to 34 Å are retrieved more efficiently than by using the "cells" command.*

>show 19

ID : 711644 -- 1

RC : a=3.83 b=3.83 c=33.62 al=90.0 be=90.0 ga=90.0

CD : sys=tetragonal den=7.4(g/cc) z=0

EM: Bi3 F O11 Pb Ti3

FO : Pb Bi3 Ti3 O11 F

NM: Lead Bismuth Titanium Oxide Fluoride

AC : a=3.833 c=33.62 spgr(A)=P4/mm

RF : J. Solid State Chem.,36,349,1981

ID : B019359 -- 2

RC : a=3.84 b=3.84 c=16.64 al=96.6 be=96.6 ga=90.0

CD : sys=tetragonal spgr(CD)=I4/mmm spno=139 den=8.0(g/cc) z=2

EM: Bi4 O12 Ti3

FO : Bi4 Ti3 O12

NM: Bismuth titanium oxide ( $4 \times 3 \times 12$ )

AC : a=3.841 c=32.83 spgr(A)=I4/mmm

RF : Ark. Kemi,1,499,1949

ID : B024961 -- 3

RC : a=3.92 b=3.92 c=14.71 al=97.7 be=97.7 ga=90.0

CD : sys=tetragonal den=7.2(g/cc) z=2

EM: Bi3 Br3 Ca O4

FO : Bi3 Ca Br3 O4

NM: Bismuth calcium bromide oxide ( $3 \times 1 \times 3 \times 4$ )

AC : a=3.92 c=28.90 spgr(A)=I

RF : Z. Anorg. Allg. Chem.,248,121,1941

ID : B024960 -- 4

RC : a=3.98 b=3.98 c=14.66 al=97.8 be=97.8 ga=90.0

CD : sys=tetragonal spgr(CD)=I4/mmm spno=139 den=7.4(g/cc) z=2

EM: Bi3 Br3 O4 Sr

FO : Bi3 Sr Br3 O4

NM: Bismuth strontium bromide oxide ( $3 \times 1 \times 3 \times 4$ )

AC : a=3.976 c=28.78 spgr(A)=I4/mmm

RF : Z. Anorg. Allg. Chem.,246,115,1941

ID : B024962 -- 5

RC : a=4.04 b=4.04 c=16.23 al=97.2 be=97.2 ga=90.0

CD : sys=tetragonal spgr(CD)=I4/mmm spno=139 den=7.4(g/cc) z=2

EM: Bi3 I3 O4 Sr

FO : Bi3 Sr I3 O4

NM: Bismuth strontium iodide oxide ( $3 \times 1 \times 3 \times 4$ )

AC : a=4.043 c=31.95 spgr(A)=I4/mmm

RF : Z. Anorg. Allg. Chem.,250,173,1942

ID : 709204 -- 6

RC : a=5.40 b=5.44 c=29.05 al=90.0 be=90.0 ga=90.0

CD : sys=orthorhombic den=3.0(g/cc) z=0

EM: Bi7 Nb O21 Ti4

FO : Bi7 Ti4 Nb O21

NM: Bismuth Titanium Niobium Oxide

AC : a=5.44 b=5.40 c=29.05 spgr(A)=P63

RF : J. Less-Common Metals,48,319,1976

ID : 805515 -- 7

RC : a=5.41 b=5.45 c=16.64 al=99.4 be=90.0 ga=90.0

CD : sys=orthorhombic spgr(CD)=C2ca spno=41 den=8.0(g/cc) z=4

EM: Bi4 O12 Ti3

FO : Bi4 Ti3 O12

NM: Tetrabismuth trititanium oxide

AC : a=5.448 b=5.411 c=32.83 spgr(A)=B2cb

RF : Ferroelectrics,3,17,1971

ID : N110983 -- 8

RC : a=5.41 b=5.45 c=32.82 al=90.0 be=90.0 ga=90.0

CD : sys=orthorhombic den=8.0(g/cc) z=4

EM: Bi4 O12 Ti3

FO : Bi4 Ti3 O12

NM: Bismuth titanium oxide (4~3~12)

AC : a=5.411 b=5.449 c=32.82

RF : J. Electrochem. Soc.,116,832,1969

*As it turns out, the last entry served as a starting point for the successful refinement of the Bi-Sr-Ca-Cu-O superconductor phase.*

The set of commands a user has to remember is small and the commands themselves are self explanatory making CRYSTDAT an efficient tool in materials research.

## References

- [1] The Powder Diffraction File, JCPDS—International Centre for Diffraction Data, Swarthmore, PA 19081.
- [2] Structure Reports, Published for the International Union of Crystallography, Vols. 8-46A (1940-1981), Reidel, Dordrecht.
- [3] Crystal Data Determinative Tables, third edition, Vol. 1 (1972), Vol. 2 (1973), Vols. 3-4 (1978), Vols. 5-6 (1983). U.S. Department of Commerce, National Bureau of Standards, and the JCPDS—International Centre for Diffraction Data, Swarthmore, PA.
- [4] CRYSTDAT (1987). An Online Search and Analysis System for NBS CRYSTAL DATA. CRYSTDAT was developed jointly by CISTI's (Canada Institute for Scientific and Technical Information) CAN/SND Scientific Numeric Database Service and the NIST Crystal Data Center.
- [5] NBS CRYSTAL DATA (1987). A Magnetic Tape of Crystallographic Data Compiled by the NIST Crystal Data Center, National Institute of Standards and Technology, Gaithersburg, MD.
- [6] Mighell, A. D., and Himes, V. L., Compound Identification and Characterization Using Lattice-Formula Matching Techniques, Acta Cryst. A42, 101 (1986).